

### CYDAL+ ALCOHOL WIPES (180'S) (CY 2003)

Each **CYDAL+ WIPE** is saturated with USP grade 70% v/v Isopropyl Alcohol and low-linting non-woven fabric. It is fast and effective against a broad spectrum of microorganisms. This is ideal for all healthcare professionals who require an effective solution for cleaning and disinfecting surfaces. Each pre-moistened wipe measures 150mm x 165mm.

#### INTENDED USES:

To clean and disinfect surfaces of non-invasive medical devices within healthcare and environment.

#### IDEAL FOR:

- Disinfecting surfaces on non-invasive medical devices and equipment.
- Use where there is high risk of infection.
- All healthcare professionals.

#### INSTRUCTIONS OF USING:

- Gently remove the cap by lifting upwards.
- Remove internal seal and feed wipe through the cap.
- Re-place the cap.
- Pull out wipes slowly from container and tear off at perforations.
- Tightly close flap between uses to retain moisture.

#### PRIMARY PACKAGING

- 180 pieces/bottle · 12 bottles/carton

#### WARNINGS:

- For professional use.
- Keep out of reach of children.
- Store below 30°C to avoid alcohol evaporation.

#### STORAGE CONDITIONS:

The product shall be stored in dry, cool and clean location. Avoid from sunlight and extreme temperature.

**MDA REGISTRATION NO.: GB6832522-99209**



**NUR JAYA MEDIC SDN. BHD. (856525-A)**

SALES AND MARKETING OFFICE:  
No. 39, Jalan PJU 1A/16, Taman Perindustrian Jaya,  
47200 Petaling Jaya, Selangor Darul Ehsan, Malaysia.  
Tel: +603-7842 8641 | Fax: +603-7842 8642

MANUFACTURING SITE (PLANT):  
No. 2 & 2A, Lorong Industri Ringan Permatang Tinggi 12,  
Taman Industri Ringan Permatang Tinggi,  
14100 Simpang Ampat, Pulau Pinang, Malaysia.  
Tel: +604-502 2474

*Quality Guaranteed. Create, Enhance & Sustain.*

✉ sales@nurjayamedic.com  
🌐 www.nurjayamedic.com

## INDEPENDENTLY TESTED AGAINST

Test Organism	Test Protocol	Log Reduction	% R	Contact Time
<i>E. hirae</i> ATCC 10541	EN 1276:2009 (E)	>5.31 ± 0.05	99.999	30 Secs
<i>E. coli</i> ATCC 10536	EN 1276:2009 (E)	>5.53 ± 0.05	99.999	30 Secs
<i>P. aeruginosa</i> ATCC 15442	EN 1276:2009 (E)	>5.37 ± 0.05	99.999	30 Secs
<i>S. aureus</i> ATCC 6538	EN 1276:2009 (E)	>5.19 ± 0.05	99.999	30 Secs
<i>Staphylococcus aureus</i> (Methicillin-resistant) ATCC 33592	EN 1276:2009 (E)	>5.40 ± 0.05	99.999	30 Secs
<i>Candida albicans</i> ATCC 10231	EN 1650:2019 (E)	>4.37 ± 0.06	99.99	30 Secs
<i>Enterococcus hirae</i> ATCC 10541	EN 13727:2012+A2:2015	>5.54 ± 0.15	99.999	30 Secs
<i>Pseudomonas aeruginosa</i> ATCC 15442	EN 13727:2012+A2:2015	>5.55 ± 0.15	99.999	30 Secs
<i>Staphylococcus aureus</i> ATCC 6538	EN 13727:2012+A2:2015	>5.49 ± 0.15	99.999	30 Secs
<i>Candida albicans</i> ATCC 10231	EN 13624:2013	>4.06 ± 0.14	99.99	30 Secs
<i>Mycobacterium Avium</i> ATCC 15679 <i>Mycobacterium terrae</i> ATCC 15755	EN 14348:2005	lg R: 6.29 lg R: >6.53	99.9999 99.9999	30 Secs
<i>Modified vaccinia virus (MVA)</i> ATCC VR-1508 <i>*According to EN 14476, a surface disinfectant is considered to possess a virucidal activity against enveloped viruses such as human coronavirus if it demonstrates a reduction of ≥4.00 log against vaccinia virus. Therefore, Nur Care Alcohol Wipes has demonstrated a virucidal activity against enveloped viruses such as human coronavirus conforming to EN 14476:2013+A2:2019 under the aforementioned conditions.</i>	EN 14476:2013+A2:2019	4.00 ± 0.00	99.99	30 Secs