

### CYDAL+ ALCOHOL FREE DISINFECTANT SPRAY (CY 2011)



**CYDAL+ ALCOHOL-FREE DISINFECTANT SPRAY** is an alcohol-free formula which comprises a synergistic mix of quaternary ammonium compound and polymeric biguanide to deliver a powerful cleaning and disinfecting solution. Ideal for all healthcare professionals who require an effective solution for cleaning and disinfecting surfaces.

#### IDEAL FOR:

- Disinfecting surfaces on non-invasive medical devices and equipment.
- Use where there is high risk of infection.
- All healthcare professionals.

#### INTENDED USES:

To clean and disinfect surfaces of non-invasive medical devices within healthcare and environment.

#### WARNING:

- For professional use.
- For external use only.
- Keep out of reach of children.

**STORAGE CONDITION:** The product shall be stored in dry, cool and clean location. Avoid from sunlight and extreme temperature.

**DISPOSAL:** Do not reuse or refill this container except as described in the directions of for use.

#### FIRST-AID MEASURES:

- Skin Contact: Wash immediately with plenty of soap and water.
- Eye Contact: Rinse continuously with water for several minutes. If irritation persists, obtain or seek medical advice.
- Ingestion: If swallowed, rinse mouth with water.
- Inhalation: Remove casualty from exposure from exposing, ensuring one's safety while doing so. Consult a doctor.

**PRIMARY PACKAGING:** 24 bottles/carton

**MDA REGISTRATION NO.:** GB7156022-97804

**NUR JAYA MEDIC SDN. BHD. (856525-A)**

SALES AND MARKETING OFFICE:  
No. 39, Jalan PJU 1A/16, Taman Perindustrian Jaya,  
47200 Petaling Jaya, Selangor Darul Ehsan, Malaysia.  
Tel: +603-7842 8641 | Fax: +603-7842 8642

MANUFACTURING SITE (PLANT):  
No. 2 & 2A, Lorong Industri Ringan Permatang Tinggi 12,  
Taman Industri Ringan Permatang Tinggi,  
14100 Simpang Ampat, Pulau Pinang, Malaysia.  
Tel: +604-502 2474

 sales@nurjayamedic.com  
 www.nurjayamedic.com

*Quality Guaranteed. Create, Enhance & Sustain.*