

CYDAL+ ALCOHOL SWAB (200'S) STERILE (CY 2001)

70% ISOPROPYL ALCOHOL

- With non woven swab
- Fabric Size:
Folded: 30mm x 30mm
Unfolded: 60mm x 30mm
- Sachet Size: 50mm x 50mm
- Thickness: 50 gsm
- 2 swabs per strip
- 200 strips per inner box (200/box)
- 50 inner boxes per outer carton

CHARACTERISTICS:

- Soft and sanitary
- Packaging is easy to tear and use
- Contains 0.5ml alcohol to maintain wetness of the pad

QUALITY CONTROL OF PRODUCT:

- The product is random checked by in-process controls and final inspected on appearance, measurements, packing etc.
- Quality and packing of the alcohol as per the general requirements for medical use.
- The Alcohol Swabs of our company are manufactured with pure non-woven fabric, with clean appearance, no fluorescence, no foreign fiber and foreign substrate.

INTENDED USE:

The Alcohol Swabs are intended for preparation of the skin by disinfecting the injection site prior to injection or venipuncture. For professional and hospital use. For external use only.

STORAGE CONDITION:

The product shall be stored in dry, cool and clean location.

INSTRUCTION FOR USE:

- Take pad out of package
- Wipe injection site vigorously
- Discard after single use
- Do not use if package is opened or damaged

MDA REGISTRATION NO.: GA6683422-93414



NUR JAYA MEDIC SDN. BHD. (856525-A)

SALES AND MARKETING OFFICE:
No. 39, Jalan PJU 1A/16, Taman Perindustrian Jaya,
47200 Petaling Jaya, Selangor Darul Ehsan, Malaysia.
Tel: +603-7842 8641 | Fax: +603-7842 8642

MANUFACTURING SITE (PLANT):
No. 2 & 2A, Lorong Industri Ringan Permatang Tinggi 12,
Taman Industri Ringan Permatang Tinggi,
14100 Simpang Ampat, Pulau Pinang, Malaysia.
Tel: +604-502 2474

Quality Guaranteed. Create, Enhance & Sustain.

✉ sales@nurjayamedic.com
🌐 www.nurjayamedic.com